COMPOUND ENTRY tutorial #2 New Compounds

Guided Data Capture (GDC)



This tutorial demonstrates how to add **COMPOUNDS**, which are **NOT** included in the database provided with the Guided Data Capture (GDC) software. □ The *GDC* software includes an integrated searchable database of more than 100,000 compounds. Most compounds can be selected from this database.

This tutorial describes how to add new compounds NOT included in the GDC compound database.



Search by Name, CASRN, or elemental formula

💐 Substance		
<u>H</u> elp		
CAS Registry Number:	Empirical formula (Case sensitive):	
Name:		Find
Search results:		Accept
		Cancel

Type the CASRN, Elemental Formula, and/or Name: (A) Elemental Formula

Enter the entire elemental formula or use wildcards (*) (e.g., C5H5N, (NH2)2CO, C*H10, C5H5*, *Cl3)

or (B) Name

Enter the entire name or use wildcards (*) (e.g., 2-methylpyridine, 2-methylp*, 2-*-pyridine) or (C) CASRN

Enter the CASRN; the final digit is filled automatically.

2. CLICK *Find* to complete the search

If you cannot find your compound in the list:



GDC compound database, but it is not exhaustive.



Guided Data Capture - Thermophysical and There Help	rmochemical Data		
Reference Compound Sample Mixture Image: Second state Substance Image: Second state Help Image: Second state CAS Begister Number.	Reaction Property Data Tables		
cyclohexane Name: Compound X name Search results:	Find Accept		
G	DC Solution State in the database? Yes No		
1. This dialog box will appear	2. CLICK Yes		
NOTE: The GDC program will assign a temporary CASRN, if you do not provide one.			
🗿 COMPOUN 🗟 Substance 🖉 🖓 🖉 🖓 🖓 🖓 🖓 🖓 🖓 🖓 🖓 🖓 🖓 🖓 🖓 🖓			







See the **SAMPLE Entry** tutorial to continue, or use the HELP screens provided within the GDC program.