


METADATA AND NUMERICAL DATA CAPTURE: **DENSITY** (2 component mixture)

Guided Data Capture (GDC)



This tutorial describes
METADATA AND NUMERICAL DATA CAPTURE:
for **DENSITIES (2 components)**
with the Guided Data Capture (GDC) software.

NOTE:

The tutorials proceed sequentially to ease the descriptions. **It is not necessary to enter *all* compounds before entering *all* samples, etc.**

Compounds, samples, properties, etc., can be added or modified at any time.

However, the hierarchy must be maintained (i.e., a property cannot be entered, if there is no associated sample or compound.)

The experimental data used in this example is from:

J. Chem. Eng. Data 2002, 47, 811–815

811

Density and Refractive Index at 298.15 K and Vapor–Liquid Equilibria at 101.3 kPa for Four Binary Systems of Methanol, *n*-Propanol, *n*-Butanol, or Isobutanol with *N*-Methylpiperazine

Shuda Chen,[†] Qunfang Lei, and Wenjun Fang*

Department of Chemistry, Zhejiang University, Hangzhou 310027, Zhejiang, People's Republic of China

Density and refractive index data at 298.15 K, and isobaric vapor–liquid equilibrium (VLE) measurements at 101.3 kPa were reported for four binary systems of methanol, *n*-propanol, *n*-butanol, or isobutanol with *N*-methylpiperazine. Excess molar volumes and refractive index deviations were calculated. The measurement results exhibit no azeotropes for VLE, negative values for excess molar volume, and positive deviations from ideality for the refractive index over the whole mole fraction range. Liquid-phase activity coefficients and vapor-phase fugacity coefficients were estimated taking into account the nonideal nature of the vapor and liquid phases. The VLE data were shown to be thermodynamically consistent and were correlated by the UNIQUAC liquid-phase activity coefficient model with temperature-dependent parameters.

Densities (liquid phase) for (methanol + N-methylpiperazine) at T = 298.15 K and p = 101.3 kPa

Table 3. Density ρ and Excess Molar Volume V^E for Binary Mixtures of Alkanol (1) + N-Methylpiperazine (2) at 298.15 K

x_1	ρ g·cm ⁻³	V^E cm ³ ·mol ⁻¹	x_1	ρ g·cm ⁻³	V^E cm ³ ·mol ⁻¹
Methanol (1) + NMP (2)			n-Propanol (1) + NMP (2)		
0.0498	0.899 54	-0.225	0.0483	0.898 02	-0.187
0.1050	0.899 14	-0.463	0.1012	0.895 74	-0.361
0.1491	0.898 68	-0.641	0.1395	0.894 07	-0.492
0.2052	0.897 88	-0.848	0.1465	0.893 78	-0.519
0.2503	0.897 23	-1.020	0.1950	0.891 38	-0.661
0.2996	0.896 17	-1.178	0.2439	0.888 98	-0.815
0.3488	0.894 96	-1.330	0.2936	0.886 06	-0.930
0.3989	0.893 39	-1.465	0.3473	0.882 79	-1.056
0.4492	0.891 59	-1.593	0.3914	0.879 95	-1.153
0.5006	0.889 11	-1.686	0.4476	0.875 98	-1.255
0.5618	0.885 69	-1.787	0.4913	0.872 65	-1.323
0.6001	0.882 75	-1.807	0.5452	0.868 03	-1.370
0.6497	0.878 42	-1.822	0.5934	0.863 55	-1.395
0.7052	0.872 25	-1.789	0.6424	0.858 59	-1.399
0.7521	0.865 54	-1.707	0.6922	0.853 04	-1.372
0.8001	0.856 87	-1.563	0.7420	0.846 78	-1.299
0.8490	0.845 58	-1.340	0.7941	0.839 47	-1.177
0.8978	0.831 09	-1.023	0.8425	0.831 79	-1.006
0.9490	0.811 44	-0.571	0.8954	0.822 30	-0.749
			0.9474	0.811 76	-0.419

Experimental Method:

Densities of all the samples were measured by using a vibrating tube digital densimeter, model DMA 602 (Anton Paar), thermostated with a circulating-water bath with a precision of ± 0.01 K. Refractive indices were determined

The data set considered here.

Guided Data Capture - Thermophysical and Thermochemical Data

File Edit Tools Help

Reference Compound Sample Mixture Reaction **Property** Data Entry

2002 che lei 0

- methanol
 - Sample 1 (cm,99.8x%,nc,x;99.8w%,glc)
- N-methylpiperazine
 - Sample 1 (cm,99.8x%,nc,mv,f;99.9w%,glc)
- methanol + N-methylpiperazine**

1. SELECT the *mixture* for which the property is to be entered.

2. CLICK *Property*

NOTE: The **bibliographic information**, **compound identities**, **sample descriptions**, and **mixture** were entered previously. (There are separate tutorials related to capture of this information.)

Start | http://pubs... | N\TRC ... | Guided D... | N\TRC ... | Microsoft... | Table of ... | Guided ... | 12:00 PM

Property Group selection

1. **CLICK** in the **property group** field.

Property and experimental method for methanol + N-methylpiperazine

Help

Property group:

Units:

Method of measurement:

Experimental purpose:

Comment (optional)

Cancel

Property group options:

- ☐ Critical properties
- ☐ Vapor pressure; Boiling temperature; and Azeotropic T & P
- ☐ Phase transition properties
- ☐ Composition at phase equilibrium
- ☐ Activity; Fugacity; and Osmotic properties
- ☒ Volumetric properties
- ☐ Heat capacity and derived properties
- ☐ Excess; partial; and apparent energetic properties

2. **SELECT volumetric properties** from the menu.

Property selection

1. **CLICK** in the **property** field.

Property and experimental method for methanol + N-methylpiperazine

Help

Property group: Volumetric properties

Property: **Specific density**

Units: Molar density

Method of measurement: 2nd virial coefficient

Experimental purpose:

Comment (optional)

Cancel

2. **SELECT specific density** from the menu.

Units selection

1. **CLICK** in the **units** field.

Property and experimental method for methanol + N-methylpiperazine

Help

Property group: Volumetric properties

Property: Specific density

Units: **kg/m³**
g/cm³
lb/ft³
ALL OTHER UNITS

Method of measurement:

Experimental purpose:

Comment (optional)

Table of values

Single value

Cancel

2. **SELECT** the **units** from the menu.

Method selection

1. CLICK in the **method of measurement** field.

Property and experimental method for methanol + N

Help

Property group: Volumetric properties

Property: Specific density

Units: g/cm³

Method of measurement:

Experimental purpose:

- Pycnometric method
- Buoyancy method
- Vibrating tube method
- Anchoring PVT measurement
- Other PVT measurement
- Burnett expansion technique
- Constant-volume piezometry
- Other experimental method (please, describe in "Comments")

Comment (optional)

Table of values

Single value

Cancel

2. SELECT a **method** from the menu, which best describes your experiment

Note: *Other* is an option. A one sentence description or a reference is often adequate.

Method detail selection

Experiment details

Select the statements, which are true for the reported measurement

- More than two calibration points used

Accept

Note: For a few methods, additional details are requested. **SELECT** those statements that apply, and **CLICK** *Accept*.

Experimental purpose selection

1. CLICK in the experimental purpose field.

Property: Specific density

Units: g/cm3

You are entering the data:

☒ In original units (as in the source) ☐ In system units (converted)

Method of measurement: Vibrating tube method

Experimental purpose:

- Principal objective of the work
- Principal objective of the work
- Secondary purpose (by-product of other objective)
- Determined for identification of a synthesized compound

Comment (optional)

2. SELECT the purpose from the menu.

Form is complete...

Property and experimental method for methanol + N-methylpiperazine

Help

Property group: Volumetric properties

Property: Specific density

Units: g/cm3

Method of measurement: Vibrating tube method Details...

Experimental purpose: Principal objective of the work


Comment (optional)

1. CLICK
Table of values

Table of values

Single value

Cancel



Specification of phases and constraints

Specific density as function of 1 variable(s)

Mixture: methanol + N-methylpiperazine

Phases in equilibrium: 1 Constraints: 0 Independent variables: Property set # 1 Sample # 1 Sample # 1 Precision of the Property Value(s) g/cm³ %

Phase of the Property Value(s)

1. SELECT # of **phases** from the pulldown list:
(one here - liquid)

2. SELECT # of **constraints** from the pulldown list:
(two here \tilde{n} T and p)

Definition of Measurement Results (Absolute vs Relative)

Data presentation
Experimental values

Comments (Optional):

Property and method Numerical Data Cancel

Specific density (kg/m³) as function of 1 variable(s)

Mixture: methanol + N-methylpiperazine

Phases in equilibrium: 1 Constraints: 2 Independent variables: 1 Property set #: 1 Sample #: 1 Sample #: 1

Phase of the Property Value(s): Precision of the Property Value(s): kg/m³ %

Definition of Measurement Results (Absolute vs Relative)

Data presentation: Experimental values

Comments (Optional):

Property and method Numerical Data Cancel

NOTE: # of variables is calculated automatically based upon the Gibbs phase rule.

NOTE: If there is more than one sample for the components, the sample # can be specified here. *(This is rare.)*

Phase of the property value selection

Specific density (g/cm³) as function of 1 variable(s)

Mixture: methanol + N-methylpiperazine

Phases in equilibrium: 1 Constraints: 2 Independent variables: 1 Property set # 1 Sample # 1 Sample # 1

Phase of the Property Value(s):

- Glass
- Glass of pure methanol
- Glass of pure N-methylpiperazine
- Liquid
- Liquid of pure methanol
- Liquid of pure N-methylpiperazine
- Liquid mixture 1
- Liquid mixture 2

Precision of the Property Value(s): g/cm³ %

Data presentation: Experimental values

Comments (Optional):

Property and method Numerical Data Cancel

1. CLICK in the Phase of the Property Values field

2. SELECT the Phase associated with the property values from the menu; (*liquid* here).

Specific density (kg/m³) as function of 1 variable(s)

Mixture:

Phases in equilibrium: Constraints: Independent variables: Property set # Sample # Sample #

Phase of the Property Value(s) Precision of the Property Value(s) ☒ kg/m³ ☐ %

Constraint 1 (Fixed value of) of Units: Uncertainty:

Constraint 2 (Fixed value of) of Units: Uncertainty:

Independent variable 1 of Units: Uncertainty: ☐ %

Definition of Measurement Results (Absolute vs Relative)

Data presentation

Comments (Optional):

NOTE: Fields for specification of *constraints, variables, units, and precisions*, appear automatically

Variable and constraint identification

Specific density (g/cm³) as function of 1 variable(s)

Mixture:

Phases in equilibrium: Constraints: Independent variables: Property set # Sample # Sample #

Phase of the Property Value(s) Precision of the Property Value(s) ☒ g/cm³ ☐ %

Constraint 1 (Fixed value of) of Value: Units: Uncertainty: ☐ %

Constraint 2 (Fixed value of) of Value: Units: Uncertainty: ☐ %

Independent variable 1 of Units: Uncertainty: ☐ %

Data presentation

Comments (Optional):

Identify constraints and variables from the menus provided.

Entry of values for constraints

Specific density (g/cm³) as function of 1 variable(s)

Mixture:

Phases in equilibrium: Constraints: Independent variables: Property set #: Sample #: Sample #:

Phase of the Property Value(s):

Precision of the Property Value(s): ☒ g/cm³ ☐ %

Constraint 1 (Fixed value of):
 of Value: Units: Uncertainty: %

Constraint 2 (Fixed value of):
 of Value: Units: Uncertainty: %

Independent variable 1:
 of

Definition of Measurement Results (Absolute vs Relative):

Data presentation:

Comments (Optional):

TYPE values for the constraints and SELECT units from the menus.

Capture of precisions, if known

The screenshot shows a software window titled "Specific density (g/cm3) as function of 1 variable(s)". The interface includes several input fields and dropdown menus. Red boxes and arrows highlight precision-related fields:

- A red box at the top right contains the label "Precision of the Property Value(s)" and a radio button selected for "g/cm3".
- Three red boxes on the right side, each with a red arrow pointing to it from a central yellow box, contain precision input fields for "Temperature", "Pressure", and "Mole fraction of methanol". The "Temperature" field contains the value "0.01".

The main interface fields include:

- Mixture: methanol + N-methylpiperazine
- Phases in equilibrium: 1
- Constraints: 2
- Independent variables: 1
- Property set #: 1
- Sample #: 1
- Phase of the Property Value(s): Liquid
- Constraint 1 (Fixed value of): Temperature of Liquid, Value: 298.15, Units: K, Uncertainty: 0.01 %
- Constraint 2 (Fixed value of): Pressure of Liquid, Value: 1, Units: bar, Uncertainty: %
- Independent variable 1: Mole fraction of methanol of Liquid, Units: Dimensionless, Uncertainty: %
- Definition of Measurement Results (Absolute vs Relative):
- Data presentation: Experimental values
- Comments (Optional):
- Buttons: Property and method, Numerical Data, Cancel

TYPE values for the precisions of the properties, variables, and constraints, if known.

NOTE: Only the precision of temperature was provided by the authors in this example.

Define results: Absolute vs Relative

Specific density (g/cm³) as function of 1 variable(s)

Mixture: methanol + N-methylpiperazine

Phases in equilibrium: 1 Constraints: 2 Independent variables: 1 Property set #: 1 Sample #: 1 Sample #: 1

Phase of the Property Value(s): Liquid Precision of the Property Value(s): g/cm³

Constraint 1 (Fixed value of): Temperature of Liquid Value: 298.15 Units: K Uncertainty: 0.01 %

Constraint 2 (Fixed value of): Pressure of Liquid Value: 1 Units: bar Uncertainty: %

Independent variable 1: Mole fraction of methanol of Liquid

Definition of Measurement Results (Absolute vs Relative)

- Direct value
- Rel. to ref. phase with same composition at FIXED T and P
- Rel. to ref. phase with SAME composition, T and P
- Rel. to mixture in EQUIL. with primary phase at same T and P
- Rel. to pure COMPONENTS in same proportion at same T and P
- Rel. to pure SOLVENT at T of same phase EQUILIBRIUM
- Rel. to pure SOLVENT at same T and P
- Rel. to pure SOLUTE at same T and P

Comments (Optional):

Property and method Numerical Data Cancel

SELECT *Direct values* from the menu.

NOTE: Other options are typically for data reported as relative values (e.g., relative to the density of the pure solvent). These are not common.

Form is complete...

Specific density (g/cm³) as function of 1 variable(s)

Mixture:

Phases in equilibrium: Constraints: Independent variables: Property set # Sample # Sample #

Phase of the Property Value(s) Precision of the Property Value(s) ☒ g/cm³ ☐ %

Constraint 1 (Fixed value of) of Value: Units: Uncertainty: ☐ %


Constraint 2 (Fixed value of) of Value: Units: Uncertainty: ☐ %

Independent variable 1 of Units: Uncertainty: ☐ %

Definition of Measurement Results (Absolute vs Relative)

Data presentation

Comments (Optional):



Specific density (g/cm³) as function of 1 variable(s)

File Edit Action Help

	Var 1	Property
1	0.0498	0.89954
2	0.1050	0.89914
3	0.1491	0.89868
4	0.2052	0.89788
5	0.2503	0.89723
6	0.2996	0.89617
7	0.3488	0.89496
8	0.3989	0.89339
9	0.4492	0.89159
10	0.5006	0.88911
11	0.5618	0.88569
12	0.6001	0.88275
13	0.6497	0.87842
14	0.7052	0.87225
15	0.7521	0.86554
16	0.8001	0.85687
17	0.8490	0.84558
18	0.8978	0.83109
19	0.9490	0.81144



Table 3. Density ρ and Excess Molar Volume V^E for Binary Mixtures of Alkanol (1) + *N*-Methylpiperazine (2) at 298.15 K

x_1	ρ	V^E	x_1	ρ	V^E
	g·cm ⁻³	cm ³ ·mol ⁻¹		g·cm ⁻³	cm ³ ·mol ⁻¹
Methanol (1) + NMP (2)			<i>n</i> -Propanol (1) + NMP (2)		
0.0498	0.899 54	-0.225	0.0483	0.898 02	-0.187
0.1050	0.899 14	-0.463	0.1012	0.895 74	-0.361
0.1491	0.898 68	-0.641	0.1395	0.894 07	-0.492
0.2052	0.897 88	-0.848	0.1465	0.893 78	-0.519
0.2503	0.897 23	-1.020	0.1950	0.891 38	-0.661
0.2996	0.896 17	-1.178	0.2439	0.888 98	-0.815
0.3488	0.894 96	-1.330	0.2936	0.886 06	-0.930
0.3989	0.893 39	-1.465	0.3473	0.882 79	-1.056
0.4492	0.891 59	-1.593	0.3914	0.879 95	-1.153
0.5006	0.889 11	-1.686	0.4476	0.875 98	-1.255
0.5618	0.885 69	-1.787	0.4913	0.872 65	-1.323
0.6001	0.882 75	-1.807	0.5452	0.868 03	-1.370
0.6497	0.878 42	-1.822	0.5934	0.863 55	-1.395
0.7052	0.872 25	-1.789	0.6424	0.858 59	-1.399
0.7521	0.865 54	-1.707	0.6922	0.853 04	-1.372
0.8001	0.856 87	-1.563	0.7420	0.846 78	-1.299
0.8490	0.845 58	-1.340	0.7941	0.839 47	-1.177
0.8978	0.831 09	-1.023	0.8425	0.831 79	-1.006
0.9490	0.811 44	-0.571	0.8954	0.822 30	-0.749
			0.9474	0.811 76	-0.419

The data set considered here.

Type or (preferably) Paste
from WORD, EXCEL, etc.

Graphical check for typographical errors


The screenshot shows a software window titled "Specific density (g/cm3) as function of 1 variable(s)". The window contains a table with 19 rows of data. The table has three columns: "Var 1", "Property", and an empty column. The data is as follows:

	Var 1	Property	
1	0.0498	0.89954	
2	0.1050	0.89914	
3	0.1491	0.89868	
4	0.2052	0.89788	
5	0.2503	0.89723	
6	0.2996	0.89617	
7	0.3488	0.89496	
8	0.3989	0.89339	
9	0.4492	0.89159	
10	0.5006	0.88911	
11	0.5618	0.88569	
12	0.6001	0.88275	
13	0.6497	0.87842	
14	0.7052	0.87225	
15	0.7521	0.86554	
16	0.8001	0.85687	
17	0.8490	0.84558	
18	0.8978	0.83109	
19	0.9490	0.81144	

At the bottom of the window, there are four buttons: "Clear the Table", "View plot", "Accept", and "Cancel". The "View plot" button is highlighted with a red rectangular box. A red curved arrow points from a yellow callout box to the "View plot" button. The callout box contains the text:


CLICK
View plot

Inspect plot




**CLICK *OK*,
if satisfied**

NOTE: The plots have several useful features. See the *HELP* screen on the plot.



Final acceptance



Guided Data Capture - Thermophysical and Thermochemical Data

File Edit Tools Help

Reference Compound Sample Mixture Reaction Property Data Tables

- 2002 che lei 0
 - methanol
 - Sample 1 (cm,99.8x%,nc;x,99.8w%,glc)
 - N-methylpiperazine
 - Sample 1 (cm,99.8x%,nc;x,99.8w%,glc)
 - methanol + N-methylpiperazine
 - ^1: VDN (Set 1), B Method:VIBTUB dT=0.01

The new dataset appears in the tree under the appropriate mixture.

Guided Data Capture - Thermophysical and Thermochemical Data

File Edit Tools Help

Reference Compound Sample Mixture Reaction

- 2002 che lei 0
 - methanol
 - Sample 1 (cm,99.8x%,nc;x;99.8w%,glc)
 - N-methylpiperazine
 - Sample 1 (cm,99.8x%,nc;mv,rf;99.9w%,glc)
 - methanol + N-methylpiperazine
 - ^1: VDN (Set 1), B Method:VIBTUB dT=0.01

NOTE: DOUBLE CLICKING on the dataset allows editing of all entered information.

END

**Continue with other compounds,
samples, properties, reactions, etc...**

***or* save your file and exit the program.**